

SenSura® Click Magnum

Med extra hög kapacitet



SenSura® Click Magnum

Med extra hög kapacitet för höga flöden



Med extra hög kapacitet

- Större volym
- Förbättrat och större utlopp
- Förbättrad diskretion
- Ökad komfort med nytt non-woven tyg
- Filter som renar odör/gaser

Produkten passar offererade SenSura Click basplattor

Kombinationen av stomivård och SenSura 2-lagershäfta har visat sig¹:

- Minska läckageproblem
- Förbättra hudkvalitén
- Förbättra livskvalitén



SenSura Magnum påse kan med fördel användas ihop med Coloplast sängpåse (träs utanpå utloppet vilket ger maximal tömning).



Sortiment

SenSura Click Magnum tömbar med backventil	Ringstorlek	Apoteksvarurn	Artikelnr	Styck/fp	Volym
	40 mm	732209	190400	10	1160 ml
	50 mm	732210	190410	10	1120 ml
	60 mm	732211	190420	10	1080 ml

1: ECET 2009 poster, "Factors affecting peristomal skin conditions – results from the Danish sub-study", 2009. WCET/WOCN 2010 Poster: Peristomal skin conditions and stoma appliances – clinical study results, O Samai et al). WCET/WOCN 2010 Poster: Factors impacting the quality of life of people with an ostomy – The DialogueStudy in the United States, J Stoia-Davis et al).